

Materials Provided:

| Name | Qty |
|--|-----|
| NanoEX™ main system | 1 |
| Power cable | 1 |
| Power adaptor | 1 |
| Serial port line | 1 |
| USB drive (For software and protocol update) | 1 |
| Sample and wash column adapters | 2 |
| Waste bottle and tubing | 1 |
| Waste bottle base | 1 |
| Accessory kit | 1 |

Materials Required but Not Provided:


| Name |
|----------------------|
| 15ml conical tube |
| Serological pipettes |

Product Specification

| Name | Spec |
|-------------------------|------------------------|
| Dimension (main system) | 22"W x 10.5"D x 14.5"H |
| Power input | 100-240V AC, 50/60Hz |

Installation Procedure

1. Connect the power cable to the power adaptor, and then connect the power adaptor to the main system.
2. Connect the power cable to a power socket then power on the system (power switch is located at the lower left-side panel).



3. Place the waste bottle on the waste bottle base and place on the right side of the main system. Securely tighten the bottle cap.
4. Connect the waste bottle tubing to the bottom of the waste connector on the main system.

